

Obrazec za popis raziskovalne opreme/ Research equipment inventory form

naziv opreme/ name of equipment	Inkubator s stresalnikom (MIR-153+MIR-S100 TABLE) MIR-153 INCUBATOR+MIR-S100 TABLE SHAKER
letnik izdelave / year of manufacture	2006
tehnični opis opreme in uporaba pri raziskavah/ technical description of equipment and use in research	<p>Inkubator je primeren za širok spekter aplikacij (inkubacija mikrobnih kultur, poskusi določanja tvorbe biofilma, protimikrobne aktivnosti) z možnostjo nastavitve temperature v razponu -10°C to 50°C. Dodana vrednost je še trostopenjski nastavljen program in stresalnik s stresalno površino 380x250 mm in razponom stresa 50-500 rpm.</p> <p>Cooled series incubator is suitable for a wide range of applications by providing a temperature range of -10°C to 50°C. In addition, MIR-153 can be applied to a wide variety of experimentation patterns with the aid of a 3-step microcomputer program. The added value is also removable shaker with a shaking surface of 380x250 mm and a range of 50-500 rpm.</p>
Pri katerem projektu je oprema bila ali se trenutno uporablja/ In which project was the equipment used or is currently in use	raziskovalno delo v sklopu ARRS projektov ali tujih partnerjev research work financed by Slovenian Research agency or partners from abroad
Lokacija opreme / Location of equipment	UP-FVZ bakteriološki laboratorij UP-FVZ bacteriology lab
Kdo opremo uporablja (raziskovalec ali oddelek)/ Who of using equipment (researcher or department)	raziskovalci/zaposleni na UP-FVZ UP-FVZ researchers/employees
Fotografije opreme/ equipment photos	



**Kontaktna oseba za opremo/
Contact person for the
equipment**

Katja Bezek

**Dodatne opombe/ Additional
notes**

<https://static1.squarespace.com/static/5876205e2e69cfaf7c9fe5a4/t/5a576cb6e2c483c38b36ee58/1515678903878/Sanyo+MIR+153-553+Instruction+Manual.pdf>

<https://www.yagev.com/media/k2/attachments/MIR-S100.pdf>