

Obrazec za popis raziskovalne opreme/ Research equipment inventory form

naziv opreme/ name of equipment	ActivPAL - Merilniki gibalne/telesne aktivnosti in sedentarnosti/ Sensors for assessing physical activity and sedentary behavior
letnik izdelave / year of manufacture	
tehnični opis opreme in uporaba pri raziskavah/ technical description of equipment and use in research	<p>ActivPAL je pospeškometer/senzor, ki se uporablja za vrednotenje količine sedentarnosti in gibalne aktivnosti otrok in mladostnikov, odraslih in starejših oseb. Senzor je velikosti 3x2x0,5 cm. Pritrdimo ga na stegno preiskovanca. ActivPAL senzorje se uporablja v raziskavah, katerih namen je objektivno spremljati gibanje preiskovancev skozi celoten dan oziroma teden./</p> <p>ActivPAL is an accelerometer, used to assess the amount of physical activity and sedentary behavior in different populations. The sensor's size is 3x2x0.5 cm. It is affixed to the thigh of the subject. These sensors are used in research, when the purpose is to objectively track the physical activity of the subjects throughout the day.</p>
Pri katerem projektu je oprema bila ali se trenutno uporablja/ In which project was the equipment used or is currently in use	<ol style="list-style-type: none"> 1. Promocija zdravja za slovenske policiste 2. Meritve sedentarnosti starejših odraslih 3. Validacija vprašalnika SST 4. Promocija zdravja v Luki Koper <ol style="list-style-type: none"> 1. Project "The promotion for Slovenian Police Officers" (PROZ-SIPO) 2. Measurements of sedentary behavior of the older adults 3. SST questionnaire validation 4. Health promotion in Port of Koper
Lokacija opreme / Location of equipment	Laboratorij za kineziologijo in fiziologijo/ Laboratory of Kinesiology and Physiology
Kdo opremo uporabljanja (raziskovalec ali oddelek)/ Who of using equipment (researcher or department)	Katedra za kineziologijo/ Department for Kinesiology Katedra za fizioterapijo/Department for Physiotherapy
Fotografije opreme/ equipment photos	



**Kontaktna oseba za opremo/
Contact person for the
equipment**

Nejc Šarabon

**Dodatne opombe/ Additional
notes**